

High-Speed Photography

PROBLEMS OF MOTION ARE EASILY ANALYZED WITH HIGH-SPEED CAMERAS



Above, a Hycam sequence. Framing rates range from 5 to 44,000 pictures per second.

Stretching time

In the study of rapid motion, photo-optical means have long been used to stretch the time base. This technique has made it possible to slow down action so it can be studied in detail.

An early trick in the motion-picture field was the slow-motion sequence, produced by speeding the camera up so that projection at normal speeds permitted newsreel viewers to appreciate the fine points of sporting events like diving, football, golf, etc.


In 1932, this speed range was substantially extended by the introduction of the rotating-prism type camera. In this development, the film moves through the

camera continuously and the image is transferred to the film by the action of a rotary optical element.

The Hycam contribution

Within the last two years, the world's most advanced rotating-prism type camera has appeared on the engineering scene and radically expanded the area of usefulness assigned to it. Going back to fundamentals, the engineers of Red Lake Laboratories Inc. used their career-long experience in high-speed photo instrumentation to refine and develop the essentially simple Hycam. In its present advanced form, Hycam exhibits the following three outstanding characteristics:

- Best picture quality by far ever produced with a rotating-prism camera (resolution—68 lines per millimeter vertically and horizontally)
- The greatest range of speeds ever available in the rotating-prism camera field (5 to 44,000 pictures per second)
- The lowest first cost, model for model, and the most economical operation, service, and maintenance.



100-ft-film-capacity Hycam—weight 13 lb 400-ft-film-capacity Hycam—weight 30 lb

more ➤

HIGH-SPEED PHOTOGRAPHY (continued)

Simplicity, ease of use, and versatility permit Hycam to be used in the laboratory by scientists, engineers, or photographers. Hycam is so designed that the following features can be added for special problems:

CINE-OSCILLOGRAPHY

In the study of mechanical equipment having electrical functions, it is often desirable to coordinate the two aspects in analysis. Replacement of the Hycam focusing eyepiece with an accessory lens permits superimpositioning of oscilloscope traces over the framing images on the film.

STREAK PHOTOGRAPHY

By replacement of the Hycam lower film gate, an


auxiliary attachment gives a slit at the film plane, adjustable from 0 to 3/16 of an inch. This permits streak photography at writing rates up to 0.08 millimeters per microsecond.

TIME-LAPSE PHOTOGRAPHY

Addition of a simple auxiliary device to Hycam permits successive exposures, individually initiated by electrical input signals at rates as low as desired. Used to speed up the action of imperceptible events.


SELF-PROJECTION

Auxiliary projection-lamp assembly converts Hycam into a projector capable of displaying any of the types of films taken by Hycam.




Oscillo Attachment

This shows the manner in which a Pentax-type lens is adapted to the camera door for cine-oscillography as described above. Weight 1 lb.



Tripods

Special tripods are available for low and high mounting of both Hycam cameras in the field. Weights 13 and 27 lb.



Timing Light Generator

Solid-state unit delivers timing pulses to neon lamp built into Hycam for putting time-base pips on the film. Operates from battery or 115-v a-c power. Weight 4 1/2 lb.




Exposure Meter

Red Lake recommends the Honeywell Pentax 1°/21° Exposure Meter, which has an angle of acceptance of only one degree and gives camera settings directly from light-level readings. Weight 18 ounces.



VISIT BOOTH 1024

Los Angeles Sports Arena
October 4-7, 1965



Power Supply

Variable autotransformer, 115-v a-c input, 0- to 320-v a-c output, controls both 100-ft and 400-ft Hycam models up to their top rated speeds. Weight 24 lb.

Red Lake Labs, Inc.

PHOTO INSTRUMENTATION EQUIPMENT

Products the User Likes



2971 CORVIN DRIVE

KIFER INDUSTRIAL PARK

SANTA CLARA, CALIFORNIA 95051

(408) 739-1698

HYCAM



1/4-20 UNC
1 HOLE


3/8-16 UNC
3 HOLES

SUB ASSY	QTY	FINAL ASSY	QTY	DRN CHK
FIRST USED ON MODEL (JOB) No:				ENG

40010

CHG
A

CHG	E. O.	BY	CHK	DATE
A	ORIG. REL.	EMW		4-17-64





EFFECTIVE FILM PLANE

WEIGHT APPROX. 13 LBS


ST'D CAMERA DIM'S SHOWN - SUBJECT TO CHANGE
WITHOUT NOTICE

INCH TOLERANCES UNLESS OTHERWISE NOTED		ITEM	PART NO	DESCRIPTION	QTY
± .050	ANGULAR ± °	FIN: —		Red Lake Laboratories SANTA CLARA, CALIF.	This document contains proprietary information and is tendered subject to the conditions that the information (A) be retained in confidence. (B) Not be reproduced or copied in whole or in part and (C) not be issued or incorporated in any product, except under an express written agreement
± .020	CONCENT. T.I.R.	H.T.: —			
± .005	MACH. SURFACES \checkmark 63/ R.	MAT'L: —			
± —	MAX. FILLET				
BREAK EDGES AND CORNERS .005— .020 DO NOT SCALE THIS DRAWING					
4-17-64	APV	OUTLINE DWG -		40010	CHG A
4-17-64	APV				
EMW	SCALE 1/3	MDL'S K2001 & K2001A HYCAM (100')			



EYE OR OSCILLO.
LENS

11 1/4




SUB ASSY	QTY	FINAL ASSY	QTY	EN
				DR
				TO
				1/8
				1/4
				1/2
				3/4
				UF
				DR
				CH
				EN

40249

CHG
A

CHG	E.O.	BY	CHK	DATE
A	ORIG. REL.	EB		5-26-65



ITEM	PART NO.	DESCRIPTION	QTY
<p>STANDARD TOLERANCES UNLESS OTHERWISE NOTED</p> <p>LED HOLES -.000 +.002 +.004 +.006 +.008 +.010</p> <p>BREAK EDGES & CORNERS .005/.010 THDS = UN CLASS 2A OR 2B DO NOT SCALE PRINT</p> <p>X ±.020 .XX ±.010 .XXX ±.005 FRACT. ± —</p> <p>ANGULAR ± CONCENT. MACH. SURF. $\frac{63}{R}$ MAX. FILL.</p> <p>T.I.R.</p>			
		Red Lake Labs, Inc.	
		PHOTO INSTRUMENTATION EQUIPMENT	
		WEIGHT, APPROX. 30 LBS	
		NOTE: STD CAMERA DIMS SHOWN- SUBJECT TO CHANGE WITHOUT NOTICE	
		SCALE 1/3	SHEET 1 OF 1
		FIRST USED ON MODEL K2004E	
		OUTLINE DRAWING	
		MODEL K2004E HYCAM(400')	
		40249	CHG A

EBred 5-26-65

ENG

APV

EMW 6-1-65

APV